

T Cell Receptor Characterization of PBMCs Using the CytoFlex

APPLICATION NOTE



Authors : Michael McPherson, Karen Carr, Si-Han Hai.

Affiliation: Beckman Coulter

IN THIS PAPER YOU WILL

Learn how to prepare PBMC and lyse red blood cells for flow cytometric staining

Learn how to stain PBMC for T cell receptor analysis by flow cytometry utilizing DuraClone

Identify TCR subpopulations using the CytoFLEX flow cytometer and DuraClone

Introduction

T cells play a pivotal role in the function and regulation of the immune system. Mature T cells express both CD3 and a T Cell Receptor (TCR). They can further be divided into two subpopulations based on their type of TCR: $\alpha\beta$ T cells (95% of the T cell repertoire) and $\gamma\delta$ T cells (5% of the T cell repertoire). TCR $\alpha\beta$ T cells are then differentiated into CD4 T cells or CD8 T cells based on their affinity recognition of antigens bound to MHC Class II or MHC Class I, respectively. The largest subset of $\gamma\delta$ T cells express the V δ 2 TCR complex. A smaller subset of $\gamma\delta$ T cells express the V δ 1 TCR complex. The study of TCR $\gamma\delta$ V δ 2 T cells and TCR $\gamma\delta$ V δ 1 T cells is important in feto-maternal¹ and allograft tolerance². The understanding of the many T cell subpopulations is important for a variety of normal and pathological studies.

DuraClone IM panels are unitized, dry format reagent cocktails that are room temperature stable. Tube formulations were based on the ONE Study as well as design input from expert flow Cytometry labs for use in clinical research studies. Here, we look at the DuraClone IM TCRs tube with normal PBMCs on the CytoFLEX.

Methods

Peripheral Blood Mononuclear Cell (PBMC) isolation

Ten milliliters of whole blood was diluted to 21 mL with 1X PBS + 3% FCS. Six mL of Ficoll (Sigma:Histopaque 1077 cat # 10771) was added to three, 15 mL conical vials. Seven mLs of diluted blood were added to the Ficoll by addition with a 10 mL pipet at a very steep angle. Diluted blood was then slowly added to over-lay on to the Ficoll. Cells were spun at 1500 RPM for 30 minutes in a swinging bucket centrifuge with the brake turned off. Supernatant was aspirate and samples were pooled and washed in PBS + 5% FCS. Cell pellets were resuspended and counted using the ViCell and diluted to a final concentration of 1×10^6 cells/100 μ L per tube for staining.

DuraClone Staining Procedure

One million PBMCs were added to a DuraClone IM TCRs tube and then vortexed for 5 seconds. After 20 minutes incubation at room temperature in the dark, the PBMCs were washed with 3 mL of PBS and centrifuged at 400g for 5 minutes. The PBMCs were resuspended in 250 μ L of PBS and acquired on the CytoFLEX.

Laser	405 nm					488 nm					638 nm		
Fluor	Krome Orange	Pacific Blue	V610	V660	V780	FITC	PE	ECD	PC5	PC7	APC	AF700	APC AF750
Marker	CD45	TCR V δ 2				TCR $\gamma\delta$	TCR $\alpha\beta$	HLA-DR		TCR V δ 1	CD4	CD8	CD3
Clone	J.33	IMMU 389				IMMU510	IP26A	Immu-357		R9.12	13B8.2	B9.11	UCHT-1

Results

Rapid TCR Subsetting with the DuraClone IM TCRs Tube

Staining of normal PBMCs with the DuraClone IM TCRs tube depicted here suggests that T cell subsets can be rapidly determined with the CytoFLEX. After approximately 30 minutes, PBMCs were accurately characterized into their appropriate compartments. Using the DuraClone IM TCRs tube is a viable alternative to laborious staining procedures.

References


- 1- Pregnancy and gamma/delta T cells: taking on the hard questions. Mincheva-Nilsson L. Reprod Biol Endocrinol. 2003 Dec 2; 1:120. Review.
- 2- Characteristics of V δ 1 (+) and V δ 2 (+) $\gamma\delta$ T cell subsets in acute liver allograft rejection. Yu X, Liu Z, Wang Y, Wang H, Zhang M, Sun Y, Su H, Jin L, Wang F, Shi M. Transpl Immunol. 2013 Dec; 29(1- 4):118-22.

Notes

The results demonstrated in this application sheet represent those generated on the Beckman Coulter CytoFLEX Flow Cytometer. As differences exist in the performance between analyzers, the author cannot guarantee a similar appearance with the use of other flow Cytometers.

Reagent Details

Reagent	Supplier	Order Details
DuraClone IM TCRs Tube	Beckman Coulter	B53340
Lymphocyte Separation Medium	Corning	25-072-C1
DPBS, 1X (Dulbecco's Phosphate-Buffered Saline)	Corning	21-031-CV



◀ **Figure 1. Rapid T Cell Receptor Subsetting with the DuraClone IM TCRs Tube**

PBMCs were gated on Lymphocytes by CD45 staining (A). CD3⁺ lymphocytes (B) were then gated based on expression of TCR $\gamma\delta$ or TCR $\alpha\beta$ (C). Of the $\gamma\delta$ T cells, expression of V δ 1 and V δ 2 was characterized (D). Of the $\alpha\beta$ T cells, CD8 and CD4 expression is characterized (E).

For Research Use Only. Not for use in diagnostic procedures.